

Agarose Anti-V5 tag antibody ab1229

★★★★★ [1 Abreviews](#) [17 References](#)

概述

产品名称	Agarose Anti-V5 tag抗体
描述	Agarose山羊多克隆抗体to V5 tag
宿主	Goat
偶联物	Agarose
经测试应用	适用于: IP
种属反应性	与反应: Species independent
免疫原	Full length native GKIPNPLLGLDST (V5) (purified): conjugated to KLH.

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

Provided as 400µl of slurry (200µl gel + 200µl buffer). 100µg of the affinity purified antibody is attached to the gel at a ratio of 500µg antibody/ml of gel or 250µg IgG/ml of 50% slurry. Affinity purified antibodies were coupled to agarose beads using a cyanogen bromide method.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.1% Sodium azide Constituents: 0.0268% PBS, 0.58% Sodium chloride
纯度	Immunogen affinity purified
克隆	多克隆

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab1229于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP	★★★★★ (1)	Use at an assay dependent concentration. Use at a concentration of 15 - 20 µl of gel slurry per 0.1 to 1 mg of protein lysate or extract.

靶标

相关性

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPNPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors