


Anti-AGBL3 antibody ab107590

2 图像

概述	
产品名称	Anti-AGBL3抗体
描述	兔多克隆抗体to AGBL3
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to amino acids 98-307 of Human AGBL3 (BC030651).
阳性对照	Lysates of: Human Fetal Liver, Human Fetal Kidney, Human Fetal Lung and Human Fetal Testis; Human Fetal Lung tissue
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Lyophilized:Reconstitute in 200ul sterile water.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG
应用	


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/2000. Predicted molecular weight: 116 kDa.
IHC-P		1/100 - 1/500.

## 靶标

功能	Metallocarboxypeptidase that may play a role in the processing of tubulin.
序列相似性	Belongs to the peptidase M14 family.
细胞定位	Cytoplasm > cytosol.

## 图片



**All lanes :** Anti-AGBL3 antibody (ab107590) at 1/1000 dilution


**Lane 1 :** Human Fetal Liver lysate

**Lane 2 :** Human Fetal Kidney lysate

**Lane 3 :** Human Fetal Lung lysate

**Lane 4 :** Human Fetal Testis lysate

**Predicted band size:** 116 kDa



Staining of Human AGBL3 in a Formalin/PFA-fixed paraffin-embedded section of Human Fetal Lung using ab107590 at a dilution of 1/100.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors