

# Anti-AKT1 (phospho S473) antibody ab8932

[60 References](#) [7 图像](#)

### 概述

产品名称	Anti-AKT1 (phospho S473)抗体
描述	兔多克隆抗体to AKT1 (phospho S473)
宿主	Rabbit
特异性	The region of AKT1 surrounding S473 has a high degree of similarity to the corresponding regions in AKT2 and AKT3 and thus may cross react with these proteins if phosphorylated on the corresponding serine residue.
经测试应用	<b>适用于:</b> ICC/IF, IHC-P, WB
种属反应性	<b>与反应:</b> Human
免疫原	Synthetic peptide corresponding to Human AKT1 aa 450-550 (C terminal) (phospho S473) conjugated to keyhole limpet haemocyanin. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.424% Potassium phosphate, 0.87% Sodium chloride
纯度	Immunogen affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using phospho peptide coupled to agarose beads followed by solid phase adsorption(s) against non-phospho peptide and non-specific peptide to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum.

克隆 多克隆  
同种型 IgG

## 应用

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab8932于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
IHC-P		1/100.
WB		Use at an assay dependent concentration. Predicted molecular weight: 56 kDa.

## 靶标

**功能** Plays a role as a key modulator of the AKT-mTOR signaling pathway controlling the tempo of the process of newborn neurons integration during adult neurogenesis, including correct neuron positioning, dendritic development and synapse formation (By similarity). General protein kinase capable of phosphorylating several known proteins. Phosphorylates TBC1D4. Signals downstream of phosphatidylinositol 3-kinase (PI(3)K) to mediate the effects of various growth factors such as platelet-derived growth factor (PDGF), epidermal growth factor (EGF), insulin and insulin-like growth factor I (IGF-I). Plays a role in glucose transport by mediating insulin-induced translocation of the GLUT4 glucose transporter to the cell surface. Mediates the antiapoptotic effects of IGF-I. Mediates insulin-stimulated protein synthesis by phosphorylating TSC2 at 'Ser-939' and 'Thr-1462', thereby activating mTORC1 signaling and leading to both phosphorylation of 4E-BP1 and in activation of RPS6KB1. Promotes glycogen synthesis by mediating the insulin-induced activation of glycogen synthase. The activated form can suppress FoxO gene transcription and promote cell cycle progression. Essential for the SPATA13-mediated regulation of cell migration and adhesion assembly and disassembly.

**组织特异性** Expressed in all human cell types so far analyzed. The Tyr-176 phosphorylated form shows a significant increase in expression in breast cancers during the progressive stages i.e. normal to hyperplasia (ADH), ductal carcinoma in situ (DCIS), invasive ductal carcinoma (IDC) and lymph node metastatic (LNMM) stages.

**疾病相关** Defects in AKT1 are a cause of susceptibility to breast cancer (BC) [MIM:114480]. A common malignancy originating from breast epithelial tissue. Breast neoplasms can be distinguished by their histologic pattern. Invasive ductal carcinoma is by far the most common type. Breast cancer is etiologically and genetically heterogeneous. Important genetic factors have been indicated by familial occurrence and bilateral involvement. Mutations at more than one locus can be involved in different families or even in the same case.  
Defects in AKT1 are associated with colorectal cancer (CRC) [MIM:114500].  
Defects in AKT1 are associated with susceptibility to ovarian cancer [MIM:604370]; also called susceptibility to familial breast-ovarian cancer type 1 (BROVCA1).

**序列相似性** Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. RAC subfamily. Contains 1 AGC-kinase C-terminal domain.

Contains 1 PH domain.  
Contains 1 protein kinase domain.


**结构域** Binding of the PH domain to the phosphatidylinositol 3-kinase alpha (PI(3)K) results in its targeting to the plasma membrane. The PH domain mediates interaction with TNK2 and Tyr-176 is also essential for this interaction.

The AGC-kinase C-terminal mediates interaction with THEM4.

**翻译后修饰** Phosphorylation on Thr-308, Ser-473 and Tyr-474 is required for full activity. Activated TNK2 phosphorylates it on Tyr-176 resulting in its binding to the anionic plasma membrane phospholipid PA. This phosphorylated form localizes to the cell membrane, where it is targeted by PDPK1 and PDPK2 for further phosphorylations on Thr-308 and Ser-473 leading to its activation. Ser-473 phosphorylation by mTORC2 favors Thr-308 phosphorylation by PDPK1. Ser-473 phosphorylation is enhanced by interaction with AGAP2 isoform 2 (PIKE-A). Ser-473 phosphorylation is enhanced in focal cortical dysplasias with Taylor-type balloon cells. Ubiquitinated; undergoes both 'Lys-48'- and 'Lys-63'-linked polyubiquitination. TRAF6-induced 'Lys-63'-linked AKT1 ubiquitination is critical for phosphorylation and activation. When ubiquitinated, it translocates to the plasma membrane, where it becomes phosphorylated. When fully phosphorylated and translocated into the nucleus, undergoes 'Lys-48'-polyubiquitination catalyzed by TTC3, leading to its degradation by the proteasome.


**细胞定位** Cytoplasm. Nucleus. Cell membrane. Nucleus after activation by integrin-linked protein kinase 1 (ILK1). Nuclear translocation is enhanced by interaction with TCL1A. Phosphorylation on Tyr-176 by TNK2 results in its localization to the cell membrane where it is targeted for further phosphorylations on Thr-308 and Ser-473 leading to its activation and the activated form translocates to the nucleus.

## 图片



ab8932 at a 1:200 dilution in Western Blot staining nuclear extract from cells infected with adenovirus expressing nuclear-targeted Akt kinase.


ab8932 at a 1:200 dilution in Western Blot staining nuclear extract from cells infected with adenovirus expressing nuclear-targeted Akt kinase.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-AKT1 (phospho S473) antibody (ab8932)

Normal human breast tissue. Akt is weakly phosphorylated in normal tissue in the breast. The phosphorylated Akt is clearly localised in the cytoplasm.

The phospho Ser 473 Akt antibody (ab8932) is used with no pretreatment at a dilution of 1:100 using Dako techmate streptavidin-biotin reagents.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-AKT1 (phospho S473) antibody (ab8932)

Akt is phosphorylated on Ser 473 in human breast tumor. The staining is much stronger than the weak basal level of phosphorylation in normal breast. The staining is cytoplasmic, as in normal tissue.


The phospho Ser 473 Akt antibody (ab8932) is used with no pretreatment at a dilution of 1:100 using Dako techmate streptavidin-biotin reagents.

This image is a magnification of the accompanying breast tumor image.




Immunocytochemistry/ Immunofluorescence - Anti-AKT1 (phospho S473) antibody (ab8932)

Confocal image of ab8932 at a 1:40 dilution staining cardiomyocytes infected with wild-type Akt in conjunction with a texas-red conjugated phalloidin to label filamentous actin in the cardiomyocytes.



Immunocytochemistry/ Immunofluorescence - Anti-AKT1 (phospho S473) antibody (ab8932)


Confocal image of ab8932 at a 1:40 dilution staining cardiomyocytes infected with wild-type Akt in conjunction with a texas-red conjugated phalloidin to label filamentous actin in the cardiomyocytes.



Immunohistochemical analysis of Human prostate cancer and melanoma cell lines, staining AKT1 (phospho S473) with ab8932.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-AKT1 (phospho S473) antibody (ab8932)

Image from Meyer S et al., PLoS One. 2012;7(6):e38222. Epub 2012 Jun 7. Fig S5.; doi:10.1371/journal.pone.0038222; June 7, 2012, PLoS ONE 7(6): e38222.



Akt is phosphorylated on Ser 473 in human breast tumor. The staining is much stronger than the weak basal level of phosphorylation in normal breast. The staining is cytoplasmic, as in normal tissue.

The phospho Ser 473 Akt antibody (ab8932) is used with no pretreatment at a dilution of 1:100 using Dako techmate streptavidin-biotin reagents.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-AKT1 (phospho S473) antibody (ab8932)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**



- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors