

Anti-alpha Lactalbumin antibody ab112972

概述

产品名称	Anti-alpha Lactalbumin抗体
描述	山羊多克隆抗体to alpha Lactalbumin
宿主	Goat
经测试应用	适用于: IHC-P, ICC, WB, ELISA
种属反应性	与反应: Cow
免疫原	The details of the immunogen for this antibody are not available.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 6.8 Preservative: 0.1% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure specificity. ab112972 was isolated by affinity chromatography using antigen coupled to agarose beads.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab112972于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 16 kDa. Colorimetric: Use 1/1000-1/10000 dilution; Chemiluminescent: Use 1/1000-1/30000 dilution.
ELISA		Use at an assay dependent concentration. Coating: Use 1/100-1/500 dilution; Primary: Use 1/1000-1/30000 dilution.

靶标

功能	Regulatory subunit of lactose synthase, changes the substrate specificity of galactosyltransferase in the mammary gland making glucose a good acceptor substrate for this enzyme. This enables LS to synthesize lactose, the major carbohydrate component of milk. In other tissues, galactosyltransferase transfers galactose onto the N-acetylglucosamine of the oligosaccharide chains in glycoproteins.
组织特异性	Mammary gland specific. Secreted in milk.
序列相似性	Belongs to the glycosyl hydrolase 22 family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors