

APC Anti-CD52 antibody [HI186], prediluted ab77025

概述

产品名称	APC Anti-CD52抗体[HI186], prediluted
描述	APC小鼠单克隆抗体[HI186] to CD52, prediluted
宿主	Mouse
偶联物	APC. Ex: 645nm, Em: 660nm
经测试应用	适用于: IHC-P, Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD52. Human tonsil
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.4 Preservative: 0.097% Sodium azide Constituents: 0.2% BSA, 0.2% PBS 0.2% BSA (high-grade protease free)
纯度	Size exclusion
纯化说明	ab77025 is purified by size exclusion chromatography and adjusted for direct use.
克隆	单克隆
克隆编号	HI186
同种型	IgG2b

The Abpromise guarantee **Abpromise™**承诺保证使用ab77025于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. ab91534 - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

靶标

功能	May play a role in carrying and orienting carbohydrate, as well as having a more specific role.
细胞定位	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors