

APC Anti-CD90 / Thy1 antibody [G7] ab25322

2 References [1 图像](#)

概述

产品名称	APC Anti-CD90 / Thy1 抗体[G7]
描述	APC大鼠单克隆抗体[G7] to CD90 / Thy1
宿主	Rat
偶联物	APC. Ex: 645nm, Em: 660nm
特异性	ab25322 recognises Thy1.1 and Thy1.2 alloantigens.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse CD90/ Thy1. Mouse T cell hybridoma C6/G8
表位	ab25322 reacts with an epitope located in the region A of CD90/Thy1.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	<p>pH: 7.3</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 16% Sucrose</p> <p>Stabilising agent.</p>
纯度	Affinity purified
克隆	单克隆
克隆编号	G7
同种型	IgG2c

轻链类型

kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab25322于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 0.2µg for 10 ⁶ cells.

靶标

功能	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

图片



Flow cytometry analysis of BALB/c mouse splenocytes stained with Rat Anti-Mouse CD90-APC ab25322 (0.01 µg/10⁶cells) and Rat Anti-Mouse CD3 epsilon-FITC **ab24948**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors