

### APC Anti-Ly6g antibody [RB6-8C5] ab25273

★★★★★ [6 Abreviews](#) [1 References](#) [1 图像](#)

#### 概述

产品名称	APC Anti-Ly6g抗体[RB6-8C5]
描述	APC大鼠单克隆抗体[RB6-8C5] to Ly6g
宿主	Rat
偶联物	APC. Ex: 645nm, Em: 660nm
经测试应用	<b>适用于:</b> Flow Cyt
种属反应性	<b>与反应:</b> Mouse
免疫原	Tissue, cells or virus corresponding to Ly6g. Mouse granulocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 16% Sucrose
纯度	Stabilising agent. Affinity purified
纯化说明	Purified from TCS
克隆	单克隆
克隆编号	RB6-8C5
同种型	IgG2b
轻链类型	kappa

## 应用

### The Abpromise guarantee

**Abpromise™**承诺保证使用ab25273于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt	★★★★★ (5)	Use a concentration of 0.25 µg/ml.

## 靶标


### 相关性

Ly6G is a GPI-anchored protein, that is also known as the myeloid differentiation antigen Gr1. The antigen is transiently expressed on monocytes in the bone marrow. The level of antigen expression in the bone marrow directly correlates with granulocyte differentiation and maturation. Ly6G is expressed predominantly on neutrophils, also in a subset of eosinophils, differentiating pre-monocytes, and plasmacytoid dendritic cells.

### 细胞定位

Attached to the membrane by a GPI-anchor.

## 图片



Flow Cytometry analysis of C57Bl/6 bone marrow cells labelling Ly6g with ab25273 at 1.125 µg (solid line) compared to a APC Rat IgG2b isotype control (dashed line) at 0.125 µg

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors