

Anti-ARHGAP18 antibody ab106553

1 References 3 图像

概述	
产品名称	Anti-ARHGAP18抗体
描述	兔多克隆抗体to ARHGAP18
宿主	Rabbit
特异性	At least two isoforms of ARHGAP18 are known to exist; ab106553 with detect both isoforms. ab106553 will not cross-react with other ARHGAP family members.
经测试应用	适用于: ICC/IF, WB, ICC
种属反应性	与反应: Mouse
免疫原	A 15 amino acid synthetic peptide near the C terminus of Human ARHGAP18 (NP_277050).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG
应用	


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 75 kDa.
ICC		Use a concentration of 2.5 µg/ml.

靶标

功能	GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.
序列相似性	Contains 1 Rho-GAP domain.

图片




**Lane 1 :** Anti-ARHGAP18 antibody (ab106553) at 1 µg/ml  
**Lane 2 :** Anti-ARHGAP18 antibody (ab106553) at 2 µg/ml

**All lanes :** 3T3 cell lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 75 kDa



ab106553 at 2.5µg/ml staining ARHGAP18 in 3T3 cells by Immunocytochemistry.

Immunocytochemistry/ Immunofluorescence - Anti-  
ARHGAP18 antibody (ab106553)

Immunofluorescence of ARHGAP18 in 3T3 cells using ab106553 at  
20 ug/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors