

# Armenian Hamster IgG monoclonal [HTK888] - Isotype Control ab18479

## 2 References

概述	
产品名称	亚美尼亚仓鼠IgG单克隆抗体[HTK888] -同型对照
经测试应用	适用于: IHC-P, IHC-Fr, IP, WB, Flow Cyt
免疫原	Chemical/ Small Molecule corresponding to Hamster. Trinitrophenol and KLH
常规说明	<p>Isotype control has been screened on a variety of resting, activated, live, and fixed mouse, rat and human tissues.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	单克隆
克隆编号	HTK888
同种型	IgG
应用	

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		
IHC-Fr		
IP		
WB		
Flow Cyt		

#### 应用说明

Flow Cyt: Use at an assay dependent dilution for surface or intracellular staining.

IHC-P: Use at an assay dependent dilution.

IHC-Fr: Use at an assay dependent dilution.

IP: Use at an assay dependent dilution.

WB: Use at an assay dependent dilution. Predicted molecular weight: n/a kDa.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors