

Anti-ASCC3 antibody ab122180

2 图像

概述

产品名称	Anti-ASCC3抗体
描述	兔多克隆抗体to ASCC3
宿主	Rabbit
经测试应用	适用于: ICC/IF, IHC-P
种属反应性	与反应: Human
免疫原	antigen sequence: KKFLNEHLQE ACTPELKPVE KTNGSFLWCE VEKYLNSTLK EMTEVPRVED LCCTLYDMLA SIKSGDELQD ELFELLGPEG LELIEKLLQN RI, corresponding to amino acids 197-288 of Human ASCC3. Run BLAST with Expasy Run BLAST with NCBI
阳性对照	Human placenta tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab122180 于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标


功能	Enhances NF-kappa-B, SRF and AP1 transactivation.
组织特异性	Ubiquitous.
序列相似性	Belongs to the helicase family. Contains 2 helicase ATP-binding domains. Contains 2 helicase C-terminal domains. Contains 3 SEC63 domains.
细胞定位	Cytoplasm.

图片



Immunofluorescent staining of Human cell line A-431 shows positivity in plasma membrane and cytoplasm. Recommended concentration of ab122180 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-ASCC3 antibody (ab122180)



ab122180 at 1/150 dilution staining ASCC3 in Paraffin-embedded Human Placenta tissue by Immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors