


Anti-Asparagine synthetase antibody [EP282Y] ab40850

重组 RabMAb

☆☆☆☆☆ 1 Abreviews 6 References 2 图像

概述

产品名称	Anti-Asparagine synthetase抗体[EP282Y]
描述	兔单克隆抗体[EP282Y] to Asparagine synthetase
宿主	Rabbit
经测试应用	适用于: WB 不适用于: Flow Cyt, ICC/IF, IHC-P or IP
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide within Human Asparagine synthetase aa 500-600 (C terminal). The exact sequence is proprietary.
阳性对照	K562 cell lysate, human breast carcinoma tissue and human muscle tissue.
常规说明	This product is a recombinant monoclonal antibody, which offers several advantages including: - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information <a href="#">see here</a> . Our RabMAb <sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a> .

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 9% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA, 50% Tissue culture supernatant
纯度	Protein A purified
克隆	单克隆

克隆编号

EP282Y

同种型

IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab40850于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★☆☆☆☆ (1)	1/1000 - 1/2000. Detects a band of approximately 60 kDa (predicted molecular weight: 61 kDa).

应用说明

Is unsuitable for Flow Cyt,ICC/IF,IHC-P or IP.

靶标


通路

Amino-acid biosynthesis; L-asparagine biosynthesis; L-asparagine from L-aspartate (L-Gln route): step 1/1.

序列相似性

Contains 1 asparagine synthetase domain.  
Contains 1 glutamine amidotransferase type-2 domain.

图片



Anti-Asparagine synthetase antibody [EP282Y] (ab40850) at 1/2000 dilution + K562 cell lysate at 10 µg

**Predicted band size:** 61 kDa  
**Observed band size:** 60 kDa

### Why choose a recombinant antibody?



**Research with confidence**  
Consistent and reproducible results



**Long-term and scalable supply**  
Recombinant technology



**Success from the first experiment**  
Confirmed specificity



**Ethical standards compliant**  
Animal-free production

Anti-Asparagine synthetase antibody [EP282Y]  
(ab40850)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors