

Anti-ATAD3B antibody ab112563

2 图像

概述

产品名称	Anti-ATAD3B抗体
描述	兔多克隆抗体to ATAD3B
宿主	Rabbit
特异性	ab112563 is predicted to not cross-react with other ATAD protein family members.
经测试应用	适用于: ICC, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human ATAD3B (C terminal). Synthetic peptide (14 amino acids) near the C terminal of Human ATAD3B (NP_114127). Database link: Q5T9A4
阳性对照	Human kidney tissue lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab112563 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab 评论	说明
ICC		Use a concentration of 5 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 73 kDa.

靶标

相关性

ATAD3B belongs to the AAA ATPase family. ATPases are a class of enzymes that catalyze the decomposition of adenosine triphosphate (ATP) into adenosine diphosphate (ADP) and a free phosphate ion. This dephosphorylation reaction releases energy, which the enzyme (in most cases) harnesses to drive other chemical reactions that would not otherwise occur. There are three named isoforms.

图片



Lane 1 : Anti-ATAD3B antibody (ab112563) at 1 µg/ml

Lane 2 : Anti-ATAD3B antibody (ab112563) at 2 µg/ml

All lanes : Human kidney tissue lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 73 kDa



Immunohistochemistry of ATAD3B in human kidney tissue with ab112563 at 5 µg/mL.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors