

Anti-ATG101 antibody ab105387

[4 References](#) [1 图像](#)

概述

产品名称	Anti-ATG101抗体
描述	兔多克隆抗体to ATG101
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	16 amino acid peptide near the center of human ATG101 (GenBank accession no. NP_001092143).
阳性对照	Human liver tissue lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab105387于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 25 kDa.

靶标

功能	Autophagy factor required for autophagosome formation.
序列相似性	Belongs to the ATG101 family.
细胞定位	Cytoplasm. Preautophagosomal structure. Under starvation conditions, it is localized to punctate structures primarily representing the isolation membrane and the isolation membrane sequesters a portion of the cytoplasm resulting in autophagosome formation.

图片



Lane 1 : Anti-ATG101 antibody (ab105387) at 1 µg/ml

Lane 2 : Anti-ATG101 antibody (ab105387) at 2 µg/ml

All lanes : human liver tissue lysate

Lysates/proteins at 15 µg/ml per lane.

Predicted band size: 25 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors