


Anti-B4GALT5 antibody ab110398

2 References 1 图像

概述	
产品名称	Anti-B4GALT5抗体
描述	兔多克隆抗体to B4GALT5
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Synthetic peptide derived from internal sequence of Human B4GALT5.
阳性对照	Extracts from HeLa cells, HuvEc cells and 293 cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG
应用	


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 45 kDa.

靶标

功能	Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.
组织特异性	Ubiquitously expressed.
通路	Protein modification; protein glycosylation.
序列相似性	Belongs to the glycosyltransferase 7 family.
细胞定位	Golgi apparatus > Golgi stack membrane. Trans cisternae of Golgi stack.

图片



Western blot - Anti-B4GALT5 antibody (ab110398)

All lanes : Anti-B4GALT5 antibody (ab110398) at 1/500 dilution

Lane 1 : Extracts from HeLa cells

Lane 2 : Extracts from 293 cells

Lane 3 : Extracts from Huvec cells

Lane 4 : Extracts from Huvec cells with immunising peptide at 10 µg

Lysates/proteins at 30 µg per lane.

Predicted band size: 45 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors