

Anti-Bacillus anthracis protective antigen antibody ab21268

1 References

概述

产品名称	Anti-Bacillus anthracis protective antigen抗体
描述	兔多克隆抗体to Bacillus anthracis protective antigen
宿主	Rabbit
经测试应用	适用于: ELISA
种属反应性	与反应: Bacillus anthracis
免疫原	Synthetic peptide corresponding to 14 amino acids at the carboxy-terminus of the Anthrax protective antigen protein (Genbank accession no. AAF86457).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab21268于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		

应用说明

ELISA: Use at a concentration of 1 µg/ml.
 Note: at this concentration ab21268 will detect 10 ng of free peptide.

Not tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

Anthrax toxins are composed of three distinct proteins, a protective antigen (PA), a lethal factor (LF) and an edema factor (EF). None of these proteins are toxic by themselves and several studies indicate that the anthrax toxin has the familiar A-B enzymatic binding structure, with PA acting as the binding domain and EF and/or LF acting as the active fragments.

细胞定位

Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors