

Anti-Bacillus thuringiensis CRY1Ab Toxin antibody ab51586

4 References

概述	
产品名称	Anti-Bacillus thuringiensis CRY1Ab Toxin抗体
描述	兔多克隆抗体to Bacillus thuringiensis CRY1Ab Toxin
宿主	Rabbit
经测试应用	适用于: IP, WB, ELISA
种属反应性	与反应: Bacillus thuringiensis
免疫原	Full length native protein (purified) corresponding to Bacillus thuringiensis CRY1Ab Toxin.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.30 Constituent: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	<u>Abpromise™</u> 承诺保证使用ab51586于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
IP		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 131 kDa.
ELISA		1/500 - 1/800.

靶标

相关性 CRY1A belongs to the delta endotoxin family and is one of the most studied insecticidal proteins produced by *Bacillus thuringiensis*. The crystal protein is produced during sporulation and is accumulated both as an inclusion and as part of the spore coat. The toxic segment of the protein is located in the N-terminus. Promotes colloid osmotic lysis by binding to the midgut epithelial cells of many lepidopteran larvae.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors