

Anti-beta Amyloid 1-38 antibody [EPR1878(2)] ab108253

重组 RabMAb

1 References 3 图像

概述

产品名称	Anti-beta Amyloid 1-38抗体[EPR1878(2)]
描述	兔单克隆抗体[EPR1878(2)] to beta Amyloid 1-38
宿主	Rabbit
特异性	ab108253 is specific to the cleaved peptide Abeta38.
经测试应用	适用于: WB, Sandwich ELISA
种属反应性	与反应: Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Abeta38 protein.
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% Tissue culture supernatant</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR1878(2)
同种型	IgG

应用


The Abpromise guarantee

Abpromise™ 承诺保证使用 ab108253 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/10000. Predicted molecular weight: 86 kDa.
Sandwich ELISA		Use a concentration of 1 µg/ml.

图片



All lanes : Anti-beta Amyloid 1-38 antibody [EPR1878(2)] (ab108253) at 1/1000 dilution


Lane 1 : Abeta38 protein

Lane 2 : Abeta40 protein

Lane 3 : Abeta42 protein

Lysates/proteins at 10 µg per lane.

Predicted band size: 86 kDa



ELISA of beta Amyloid 1-38 peptide at 2µg/mL. Detected with ab108253 at 0~500 pg/mL. Peroxidase Streptavidin SA-HRP at 1/20000 was used as a secondary antibody.

Washing buffer: 1X PBST

Blocking buffer: 1% BSA/PBS

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-beta Amyloid 1-38 antibody [EPR1878(2)]
(ab108253)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors