

Biotin Anti-beta 2 Defensin/BD-2 antibody ab83509

1 References

概述

| | |
|-------|---|
| 产品名称 | 生物素Anti-beta 2 Defensin/BD-2抗体 |
| 描述 | 生物素山羊多克隆抗体to beta 2 Defensin/BD-2 |
| 宿主 | Goat |
| 偶联物 | Biotin |
| 经测试应用 | 适用于: Sandwich ELISA, WB, ELISA |
| 种属反应性 | 与反应: Human |
| 免疫原 | Recombinant full length protein corresponding to Human beta 2 Defensin/BD-2. |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|--|
| 形式 | Lyophilized:Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml. |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 7.20 Constituent: PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab83509于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----------------|------|---|
| Sandwich ELISA | | Use at an assay dependent concentration. |
| WB | | Use at an assay dependent concentration. |
| ELISA | | Use at an assay dependent concentration. Can be paired for ELISA with Rabbit polyclonal to beta 2 Defensin / hBD-2 (ab9871). |

靶标

| | |
|-------|---|
| 功能 | Has antibacterial activity. |
| 组织特异性 | Expressed in the skin and respiratory tract. |
| 序列相似性 | Belongs to the beta-defensin family. LAP/TAP subfamily. |
| 细胞定位 | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors