

Biotin Anti-Borrelia burgdorferi antibody ab69158

概述

产品名称	生物素Anti-Borrelia burgdorferi抗体
描述	生物素兔多克隆抗体to Borrelia burgdorferi
宿主	Rabbit
偶联物	Biotin
特异性	ab69158 reacts with Borrelia burgdorferi in bacterial and infected tissue samples.
经测试应用	适用于: IHC, ICC/IF
种属反应性	与反应: Borrelia burgdorferi
免疫原	Tissue, cells or virus corresponding to Borrelia burgdorferi. Borrelia burgdorferi whole cell preparation.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
纯度	Affinity purified
克隆	多克隆
同种型	IgG

应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent dilution.

靶标

相关性 *Borrelia burgdorferi* is a spirochete and the cause of Lyme disease, a tick transmitted illness of humans and animals. *B. burgdorferi* may persist in humans and animals for months or years following initial infection, despite a robust humoral immune response. *B. burgdorferi* resembles other spirochetes in that it is a highly specialized, motile, two-membrane, spiral shaped bacteria which lives primarily as an extracellular pathogen. *B. burgdorferi* has an unusual genome compared with other eubacteria which includes a linear chromosome approximately one megabase in size and numerous linear and circular plasmids.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors