

Biotin Anti-CD200 / OX2 antibody [MRC OX90] ab33735

概述

产品名称	生物素Anti-CD200 / OX2抗体[MRC OX90]
描述	生物素大鼠单克隆抗体[MRC OX90] to CD200 / OX2
宿主	Rat
偶联物	Biotin
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Fusion protein corresponding to Mouse CD200/ OX2. Database link: O54901

常规说明

Spleen cells from immunised rats were fused with cells of the rat Y3 myeloma cell line.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein G purified
纯化说明	This purified antibody was prepared by affinity chromatography on Protein G from tissue culture supernatant.
Primary antibody说明	Spleen cells from immunised rats were fused with cells of the rat Y3 myeloma cell line.
克隆	单克隆
克隆编号	MRC OX90

应用

The Abpromise guarantee Abpromise™承诺保证使用ab33735于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Use 10µl of the suggested working dilution to label 10 ⁶ cells in 100µl.

靶标

功能	Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues.
序列相似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors