

Biotin Anti-CD32B + CD32A antibody [AT10] ab30355

[4 References](#) [1 图像](#)

概述

产品名称	生物素Anti-CD32B + CD32A抗体[AT10]
描述	生物素小鼠单克隆抗体[AT10] to CD32B + CD32A
宿主	Mouse
偶联物	Biotin
特异性	According to the CD32A + CD32B Cluster workshop (Leucocyte Typing V, pp823, pp828), clone AT10 reacts equally well with both Fc gamma RII a and b-expressing cells.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD32A. K562 cell line
阳性对照	Human peripheral blood monocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	AT10
骨髓瘤	NS1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab30355于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt		Use 10µl for 10 ⁶ cells. ab18434 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标

细胞定位

CD32B: Cell membrane. CD32A: Cell membrane.

图片



Flow cytometric analysis of human peripheral blood monocytes with 10µl for 10⁶ cells of anti human CD32A + CD32B Biotin conjugated (ab30355).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors