

Biotin Anti-CXCL14 antibody ab46008

1 图像

概述

产品名称	生物素Anti-CXCL14抗体
描述	生物素兔多克隆抗体to CXCL14
宿主	Rabbit
偶联物	Biotin
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Recombinant fragment
免疫原	Highly pure (>98%) recombinant hBRAK (human BRAK).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Lyophilized:Reconstitute to a concentration of 50µg/ml with sterile PBS solution containing 0.1% BSA.
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: PBS, 0.1% BSA
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab46008于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		Use a concentration of 0.25 - 1 µg/ml. Using 100 µl/well antibody solution.

靶标

功能	Potent chemoattractant for neutrophils, and weaker for dendritic cells. Not chemotactive for T-cells, B-cells, monocytes, natural killer cells or granulocytes. Does not inhibit proliferation of myeloid progenitors in colony formation assays.
组织特异性	Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Highly expressed in normal tissue without inflammatory stimuli and infrequently expressed in cancer cell lines. Weakly expressed in monocyte-derived dendritic cells. Not detected in lung or unstimulated peripheral blood lymphocytes.
序列相似性	Belongs to the intercrine alpha (chemokine CxC) family.
细胞定位	Secreted.

图片



ab46008 detecting CXCL14 by sandwich ELISA using 100 µl/well antibody solution. Samples were incubated with primary antibody ab46008 (1.0 µg/ml) and capture antibody Polyclonal Anti-Human CXCL14. Allows the detection of at least 0.2 – 0.4 ng/well of recombinant CXCL4

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors