


### Biotin Anti-Endothelin 1 antibody ab48251

#### 概述

产品名称	生物素Anti-Endothelin 1抗体
描述	生物素兔多克隆抗体to Endothelin 1
宿主	Rabbit
偶联物	Biotin
经测试应用	适用于: ELISA, RIA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Sheep, Rabbit, Pig 
免疫原	Synthetic peptide: SCSSLMNKECVYFCHLDIIW conjugated to KLH by a Cysteine residue linker, corresponding to amino acids 54-73 of Human Endothelin 1

 [Run BLAST with](#)

 [Run BLAST with](#)

#### 常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 49.73% PBS, 50% Glycerol (glycerin, glycerine), 0.25% BSA
纯度	Protein G purified
克隆	多克隆
同种型	IgG

## 应用

**The Abpromise guarantee** Abpromise™承诺保证使用ab48251于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.

## 靶标

功能	Endothelins are endothelium-derived vasoconstrictor peptides.
组织特异性	Expressed in lung, placental stem villi vessels and in cultured placental vascular smooth muscle cells.
序列相似性	Belongs to the endothelin/sarafotoxin family.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors