

Biotin Anti-IL-13 antibody [JES10-2E10] ab24955

概述

产品名称	生物素Anti-IL-13抗体[JES10-2E10]
描述	生物素大鼠单克隆抗体[JES10-2E10] to IL-13
宿主	Rat
偶联物	Biotin
经测试应用	适用于: ELISA
种属反应性	与反应: Human, Cynomolgus monkey, Rhesus monkey
免疫原	Recombinant full length protein corresponding to Human IL-13.
常规说明	<p>This is a neutralizing antibody.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
存储溶液	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Affinity purified
Primary antibody说明	This is a neutralizing antibody.
克隆	单克隆
克隆编号	JES10-2E10
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab24955于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/2000 - 1/4000.

靶标

功能	Cytokine (PubMed:8096327, PubMed:8097324). Inhibits inflammatory cytokine production (PubMed:8096327). Synergizes with IL2 in regulating interferon-gamma synthesis (PubMed:8096327). May be critical in regulating inflammatory and immune responses (PubMed:8096327, PubMed:8097324). Positively regulates IL31RA expression in macrophages.
疾病相关	Allergic rhinitis
序列相似性	Belongs to the IL-4/IL-13 family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors