

Biotin Anti-IL-1RA antibody [A71B6D11] ab25074

★★★★★ 2 Abreviews

概述

产品名称	生物素Anti-IL-1RA抗体[A71B6D11]
描述	生物素小鼠单克隆抗体[A71B6D11] to IL-1RA
宿主	Mouse
偶联物	Biotin
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human IL-1RA.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.50</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 0.268% PBS, 0.5% BSA</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	A71B6D11
骨髓瘤	NS0
同种型	IgG1

The Abpromise guarantee **Abpromise™**承诺保证使用ab25074于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA	★★★★★ (2)	Use at an assay dependent concentration.

靶标

功能	Inhibits the activity of interleukin-1 by binding to receptor IL1R1 and preventing its association with the coreceptor IL1RAP for signaling. Has no interleukin-1 like activity. Binds functional interleukin-1 receptor IL1R1 with greater affinity than decoy receptor IL1R2; however, the physiological relevance of the latter association is unsure.
组织特异性	The intracellular form of IL1RN is predominantly expressed in epithelial cells.
疾病相关	Microvascular complications of diabetes 4 Interleukin 1 receptor antagonist deficiency
细胞定位	Cytoplasm and Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors