

Biotin Anti-Melanoma Inhibitory Activity antibody ab84235

★★★★★ 1 Abreviews

概述

产品名称	生物素Anti-Melanoma Inhibitory Activity抗体
描述	生物素兔多克隆抗体to Melanoma Inhibitory Activity
宿主	Rabbit
偶联物	Biotin
经测试应用	适用于: Sandwich ELISA, Flow Cyt, WB
种属反应性	与反应: Human
免疫原	Highly pure (>98%) recombinant full length Human Melanoma Inhibitory Activity
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Lyophilized:Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: PBS
纯化说明	ab84235 was purified by affinity chromatography and then biotinylated.Filter sterilised.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab84235于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		Use at an assay dependent concentration.
Flow Cyt	★★★★★ (1)	Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

靶标

功能	Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.
组织特异性	All malignant melanoma cell lines tested and infrequently in glioma cell lines.
序列相似性	Belongs to the MIA/OTOR family. Contains 1 SH3 domain.
翻译后修饰	May possess two intramolecular disulfide bonds.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors