

# Biotin Mouse IgG3, kappa monoclonal [MG3-35] - Isotype Control ab18436

### 概述

产品名称	生物素小鼠IgG3, kappa单克隆抗体[MG3-35] -同型对照
偶联物	Biotin
特异性	The isotype control has been screened on a variety of resting, activated, live, and fixed mouse and human tissues.
经测试应用	<b>适用于:</b> Flow Cyt
免疫原	Chemical/ Small Molecule corresponding to Mouse conjugated to keyhole limpet haemocyanin. Trinitrophenol and KLH
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	单克隆
克隆编号	MG3-35
同种型	IgG3
轻链类型	kappa

### 应用

## The Abpromise guarantee

**Abpromise™**承诺保证使用ab18436于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors