

Biotin Anti-V5 tag antibody [SV5-Pk1] ab27672

概述

产品名称	生物素Anti-V5 tag抗体[SV5-Pk1]
描述	生物素小鼠单克隆抗体[SV5-Pk1] to V5 tag
宿主	Mouse
偶联物	Biotin
经测试应用	适用于: IHC-Fr, ELISA, WB
种属反应性	与反应: Species independent
免疫原	Tissue, cells or virus corresponding to V5 tag. Paramyxovirus Simian-Virus 5 (SV5)
表位	Ab27672 recognises a small epitope, termed Pk, present on the P/V proteins of the paramyxovirus, SV5.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide
纯度	Protein G purified
纯化说明	This antibody was purified from tissue culture supernatant by affinity chromatography.
克隆	单克隆
克隆编号	SV5-Pk1
骨髓瘤	Sp2/0-Ag14
同种型	IgG2a

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab27672于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
WB		1/1000 - 1/5000.

靶标

相关性 The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors