

# Anti-Blood Group AB antigen antibody [Z5H-2/Z2A] ab24223

## 2 References

### 概述

产品名称	Anti-Blood Group AB antigen抗体[Z5H-2/Z2A]
描述	小鼠单克隆抗体[Z5H-2/Z2A] to Blood Group AB antigen
宿主	Mouse
经测试应用	适用于: ELISA, IP, ICC/IF, Flow Cyt, WB
种属反应性	与反应: Human
免疫原	Antigen A B of blood cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	Z5H-2/Z2A
同种型	IgM

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab24223于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration. PubMed: 19608744
Flow Cyt		Use at an assay dependent concentration. <b>ab91545</b> - Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody.
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors