

# Anti-Bordetella pertussis toxin/ptxD antibody [56.2] ab37686

## 1 References

### 概述

产品名称	抗百日咳杆菌toxin/ptxD抗体[56.2]
描述	小鼠单克隆抗体[56.2] to百日咳杆菌toxin/ptxD
宿主	Mouse
特异性	This product reacts with different epitope than <a href="#">ab37574</a> , <a href="#">ab37684</a> and <a href="#">ab37687</a> .
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Bordetella pertussis
免疫原	<i>Bordetela pertussis</i> Toxin B-Oligomer.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C long term.
存储溶液	pH: 7.40 Constituents: 2.9% Sodium chloride, 0.0268% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	56.2
骨髓瘤	x63-Ag8.653
同种型	IgG1
轻链类型	kappa

### 应用

## The Abpromise guarantee

**Abpromise™**承诺保证使用ab37686于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Detects a band of approximately 11.7 kDa.

## 靶标

### 相关性

Pertussis toxin contains five different chains, S1-S5. They are organized into 2 functional subunits: A, composed of S1 (which is toxic) and B, containing S2, S3, S5, and two copies of S4 (B binds to the membrane receptors). Dimers of S2-S4 and S3-S4 are held together by S5

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors