

# Anti-BOTOX antibody ab20641

★★★★★ [1 Abreviews](#) [5 References](#)

### 概述

---

产品名称	Anti-BOTOX抗体
描述	兔多克隆抗体to BOTOX
宿主	Rabbit
特异性	Clostridium botulinum Type A Neurotoxin.
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Species independent
免疫原	Full length protein corresponding to BOTOX.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 2.68% PBS, 20% Glycerol (glycerin, glycerine)
纯度	Protein G purified
克隆	多克隆
同种型	IgG

### 应用

---

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab20641于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/10000.
WB	★★★★★ (1)	Use at an assay dependent dilution. Predicted molecular weight: 149 kDa.

## 靶标

**相关性** Clostridium botulinum neurotoxins are zinc-dependent proteases that block the release of acetylcholine from peripheral cholinergic nerve endings by targeting specific membrane proteins including synaptosomal-associated proteins (SNAP), vesicle-associated membrane proteins (VAMP), and syntaxins. Neurotoxin A is characteristically associated with human botulism and known to target and hydrolyse the 197-Gln-|-Arg-198 bond of SNAP-25.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors