

Anti-BTG2/PC3 antibody ab85051

★★★★☆ [1 Abreviews](#) [8 References](#) [2 图像](#)

概述

产品名称	Anti-BTG2/PC3抗体
描述	兔多克隆抗体to BTG2/PC3
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human BTG2/PC3 (C terminal). A 15 amino acid synthetic peptide from the C-terminal residues. Database link: NM_006763
阳性对照	WB: HepG2 whole cell lysate (ab7900). IHC-P: Human fetal colon tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Lyophilized:Reconstitute in 200ul sterile H2O.
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 2% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab85051 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★☆ (1)	1/200 - 1/1000. Detects a band of approximately 18 kDa (predicted molecular weight: 18 kDa).
IHC-P		1/100 - 1/500.

靶标


功能

Involved in cell cycle regulation. Could be involved in the growth arrest and differentiation of the neuronal precursors (By similarity). Anti-proliferative protein. Modulates transcription regulation mediated by ESR1. Involved in neurite outgrowth.

序列相似性


Belongs to the BTG family.

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human fetal colon tissue labelling BTG2/PC3 with ab85051 at a dilution of 1/100.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-BTG2/PC3 antibody (ab85051)



Anti-BTG2/PC3 antibody (ab85051) at 1/250 dilution + HepG2 cell lysates

Predicted band size: 18 kDa

Observed band size: 18 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors