

Anti-C4 antibody [16D2] ab11863

★★★★★ [1 Abreviews](#) [13 References](#) [1 图像](#)

概述

产品名称	Anti-C4抗体[16D2]
描述	大鼠单克隆抗体[16D2] to C4
宿主	Rat
特异性	Antibody cross reacts with C4, C4b and C4d.
经测试应用	适用于: IHC-Fr
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus. This information is proprietary to Abcam and/or its suppliers.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
纯度	Protein G purified
克隆	单克隆
克隆编号	16D2
同种型	IgG2a

应用

The Abpromise guarantee


Abpromise™ 承诺保证使用ab11863于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.

靶标

图片



Frozen section of mouse pulmonary artery of control (left) and mouse with hypoxia-induced pulmonary hypertension (right). Staining of C4 (red) with ab11863 in a dilution 1/300. Green is autofluorescence of elastic lamellae defining vascular media. Blue is cell nucleus.

Immunohistochemistry (Frozen sections) - Anti-C4 antibody [16D2] (ab11863)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors