

### Anti-C5a / C5a des Arg antibody [2952] ab11878

★★★★★ [1 Abreviews](#) [5 References](#)

#### 概述

产品名称	Anti-C5a / C5a des Arg抗体[2952]
描述	小鼠单克隆抗体[2952] to C5a / C5a des Arg
宿主	Mouse
特异性	<p>This antibody reacts with a neo-epitope (des-Arg) on C5a that is formed upon cleavage of C5 into C5a and C5b. It will detect C5a and C5a des-Arg (Not C5).</p> <p>Please notice that under given conditions it is known that C5 can expose epitopes normally only found in the cleaved activation products.</p>
经测试应用	<b>适用于:</b> ELISA, IHC-Fr, Functional Studies, WB, IHC-P
种属反应性	<b>与反应:</b> Human
免疫原	Other Immunogen Type. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Human Granulocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: PBS, 0.1% BSA
纯度	Protein G purified
纯化说明	0.2 µm filtered
克隆	单克隆
克隆编号	2952
同种型	IgG1

## 应用

### The Abpromise guarantee

**Abpromise™**承诺保证使用ab11878于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA	★★★★★ (1)	Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Functional Studies		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

## 靶标

### 相关性

C5 is involved in the activation of the lytic pathway within the complement system, which is an important factor in innate immunity.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors