

Anti-C5a-R antibody [P12/1] ab24036

1 References

概述

产品名称	Anti-C5a-R抗体[P12/1]
描述	小鼠单克隆抗体[P12/1] to C5a-R
宿主	Mouse
特异性	No P12/1 signal after pre-incubation of HMC-1 cells with C5a, meaning that P12/1 does not displace C5a.
经测试应用	适用于: Flow Cyt, WB, IHC-P, IHC-Fr
种属反应性	与反应: Human, Rhesus monkey
免疫原	Synthetic peptide corresponding to Human C5a-R aa 1-31 (N terminal). Sequence: MNSFNYTTDPYGHYDDKDTLDLNTVPVDKTSN



[Run BLAST with](#)



[Run BLAST with](#)

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	pH: 7.3 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Affinity purified
克隆	单克隆
克隆编号	P12/1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab24036于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 1 mg/ml.
WB		Use a concentration of 1 mg/ml. Predicted molecular weight: 39 kDa.
IHC-P		Use a concentration of 1 mg/ml.
IHC-Fr		Use a concentration of 1 mg/ml.

靶标

功能	Receptor for the chemotactic and inflammatory peptide anaphylatoxin C5a. This receptor stimulates chemotaxis, granule enzyme release and superoxide anion production.
序列相似性	Belongs to the G-protein coupled receptor 1 family.
翻译后修饰	Sulfation plays a critical role in the association of the receptor with C5a, but no significant role in the ability of the receptor to transduce a signal and mobilize calcium in response to a small peptide agonist.
细胞定位	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors