

Anti-Calretinin antibody [6B3] ab82198

★★★★★ 1 Abreviews

概述

产品名称	Anti-Calretinin抗体[6B3]
描述	小鼠单克隆抗体[6B3] to Calretinin
宿主	Mouse
特异性	No cross-reaction with Calbindin or Calmodulin.
经测试应用	适用于: WB
种属反应性	与反应: Rat, Human
免疫原	Recombinant Human Calretinin 22kDa
阳性对照	Cerebral cortex, pleural mesothelioma.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.30 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	6B3
同种型	IgG

The Abpromise guarantee **Abpromise™**承诺保证使用ab82198于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		

应用说明

WB: 1/40 - 1/50. Predicted molecular weight: 31 kDa.

The antibody is suitable for the detection of Calretinin in neurons and may be used as an aid in differential diagnosis of pleura-mesothelioma. This marker is frequently expressed in epitheloid mesothelioma, while occasionally it turns positive in sarcomatous mesothelioma and carcinoma. It was found to be negative in sarcoma and in 90% of adenocarcinoma.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能 Calretinin is a calcium-binding protein which is abundant in auditory neurons.

组织特异性 Brain.

序列相似性 Belongs to the calbindin family.

Contains 6 EF-hand domains.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors