

Anti-CCDC14 antibody ab122262

2 图像

概述

产品名称	Anti-CCDC14抗体
描述	兔多克隆抗体to CCDC14
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	antigen sequence: IYQALCEHV QTQMSLMNDL TSKNIPNGIP AVPCHAPSHS ESQATPHSSY GLCTSTPVWS LQRPPCPPKV HSEVQTDGNS QFASQG, corresponding to amino acids 205-290 of Human CCDC14. Run BLAST with Expasy Run BLAST with NCBI
阳性对照	Human cerebral cortex tissue; Human plasma lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee


Abpromise™ 承诺保证使用ab122262于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/1000 - 1/2500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use at an assay dependent concentration.


靶标

图片



Immunocytochemistry/ Immunofluorescence - Anti-CCDC14 antibody (ab122262)

Immunofluorescent staining of human cell line A-431 shows positivity in cytoplasm, vesicles & nucleus but excluded from the nucleoli.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CCDC14 antibody (ab122262)

ab122262, at 1/1200, staining CCDC14 in paraffin embedded Human cerebral cortex tissue by Immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors