


# Anti-CD15 antibody [28] ab20137

## 5 References

### 概述

产品名称	Anti-CD15抗体[28]
描述	小鼠单克隆抗体[28] to CD15
宿主	Mouse
特异性	Specific to CD15. 28 is useful for leukaemia typing and purification of stem cells.
经测试应用	适用于: IP, Flow Cyt, IHC-Fr
种属反应性	预测可用于: Human 
免疫原	Tissue, cells or virus corresponding to CD15. Monocytes, separated from other peripheral blood leucocytes (PBL) on fibronectin plates.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	28
骨髓瘤	Sp2/0-Ag14
同种型	IgM
轻链类型	unknown

## 应用

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab20137于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. <b>ab91545</b> - Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody.
IHC-Fr		Use at an assay dependent concentration. Note: samples should be acetone fixed.

## 靶标

**相关性** CD15 is a carbohydrate adhesion molecule (and not a protein) that mediates phagocytosis and chemotaxis. Synthesis is directed by FUT4 in lymphoid cells and mature granulocytes, and by FUT9 in promyelocytes and monocytes.

**细胞定位** Golgi Apparatus; Membrane-bound form in trans cisternae of Golgi.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors