

Anti-CD160 antibody [CNX46-3] ab91101

1 References 1 图像

概述

产品名称	Anti-CD160抗体[CNX46-3]
描述	大鼠单克隆抗体[CNX46-3] to CD160
宿主	Rat
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Recombinant Mouse CD160-Fc protein
阳性对照	Flow Cyt: Mouse (C57BL/6) splenocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	CNX46-3
同种型	IgG2a
轻链类型	kappa

The Abpromise guarantee      **Abpromise™**承诺保证使用ab91101于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. 1 µg/test

靶标

功能	Receptor showing broad specificity for both classical and non-classical MHC class I molecules.
组织特异性	Expressed in spleen, peripheral blood, and small intestine. Expression is restricted to functional NK and T cytotoxic lymphocytes.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

图片



Flow cytometry analysis of C57BL/6 splenocytes labelling CD160 with ab91101 at 0.5 µg, followed by F (ab')2 Anti-Rat IgG PE and subsequently with anti-Mouse CD3e FITC and anti-Mouse NK1-1 APC. The histogram (right) demonstrates staining of **ab911001** on NK1-1+CD3. cells (red histogram) and NK1-1+CD3+ cells (green histogram), as gated in the dot plot (left).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors