

Anti-CD1a antibody [NA1/34-HLK] ab24055

6 References

概述

产品名称	Anti-CD1a抗体[NA1/34-HLK]
描述	小鼠单克隆抗体[NA1/34-HLK] to CD1a
宿主	Mouse
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Dog, Human, Monkey
免疫原	Tissue, cells or virus corresponding to Human CD1a. Thymocytes (Human).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
纯度	Ion Exchange Chromatography
纯化说明	Purified IgG prepared by ion exchange chromatography
克隆	单克隆
克隆编号	NA1/34-HLK
骨髓瘤	NS1/1-Ag4-1
同种型	IgG2a

The Abpromise guarantee **Abpromise™**承诺保证使用ab24055于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Use 10 ul of this concentration to label 10^6 cells in 100ul. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

靶标

功能	Antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells.
组织特异性	Expressed on cortical thymocytes, epidermal Langerhans cells, dendritic cells, on certain T-cell leukemias, and in various other tissues.
序列相似性	Contains 1 Ig-like (immunoglobulin-like) domain.
细胞定位	Cell membrane. Endosome membrane. Subject to intracellular trafficking between the cell membrane and endosomes. Localizes to cell surface lipid rafts.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors