

Anti-CD27 antibody [LG.3A10] ab219779

[1 References](#) [1 图像](#)

概述

产品名称	Anti-CD27抗体[LG.3A10]
描述	亚美尼亚仓鼠单克隆抗体[LG.3A10] to CD27
宿主	Armenian hamster
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Mouse CD27. Hamster spleen cells were fused with SP2/0 murine myeloma cells. Database link: P26842
阳性对照	Human tonsil tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Proprietary Purification
克隆	单克隆
克隆编号	LG.3A10
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab219779 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 15 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

靶标

功能	Receptor for CD70/CD27L. May play a role in survival of activated T-cells. May play a role in apoptosis through association with SIVA1.
组织特异性	Found in most T-lymphocytes.
序列相似性	Contains 3 TNFR-Cys repeats.
翻译后修饰	Phosphorylated and O-glycosylated.
细胞定位	Membrane.

图片



Immunohistochemical analysis of formalin-fixed, paraffin-embedded human tonsil tissue labeling CD27 with ab219779 at 15 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD27 antibody [LG.3A10] (ab219779)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors