

Anti-CD28 antibody [JJ319] - Low endotoxin, Azide free ab185775

1 图像

概述

产品名称	Anti-CD28抗体[JJ319] - Low endotoxin, Azide free
描述	小鼠单克隆抗体[JJ319] to CD28 - Low endotoxin, Azide free
宿主	Mouse
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Rat
免疫原	The details of the immunogen for this antibody are not available.
阳性对照	Lou Rat splenocytes.
常规说明	<p>Endotoxin Level: Less than 0.001 ng/μg antibody, as determined by the LAL assay.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2
无载体	Aqueous buffer, no sodium azide.
纯度	是
克隆	Affinity purified
克隆编号	单克隆
克隆编号	JJ319

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab185775于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 1µg for 10 ⁵⁻⁸ cells. (In a final volume of 100 µL).

靶标

功能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
组织特异性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Membrane.

图片



Flow cytometric analysis of Lou Rat splenocytes with 0.5 µg Mouse IgG1 K Isotype Control (open histogram) or 0.5 µg of ab185775 (filled histogram) followed by Anti-Mouse IgG FITC. Total viable cells were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors