

Anti-CD30 antibody [MM0169-5K28] ab89194

概述

产品名称	Anti-CD30抗体[MM0169-5K28]
描述	小鼠单克隆抗体[MM0169-5K28] to CD30
宿主	Mouse
经测试应用	适用于: Flow Cyt, WB, Functional Studies
种属反应性	与反应: Human
免疫原	Human recombinant CD30 extracellular domain
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
纯化说明	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography.
克隆	单克隆
克隆编号	MM0169-5K28
同种型	IgG2

应用

The Abpromise guarantee Abpromise™ 承诺保证使用ab89194于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		
WB		
Functional Studies		

应用说明

Flow Cyt: 1/20 - 1/100.
 FuncS: CD30 activation: Use at an assay dependent dilution.
 WB: 1/500 - 1/1000. Predicted molecular weight: 64 kDa.

Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

靶标

功能

Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.

序列相似性

Contains 6 TNFR-Cys repeats.

翻译后修饰

Phosphorylated on serine and tyrosine residues.

细胞定位

Cytoplasm and Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors