


Anti-CD31 antibody ab28364

★★★★☆ [156 Abreviews](#) [2004 References](#) [3 图像](#)

概述

产品名称	Anti-CD31抗体
描述	兔多克隆抗体to CD31
宿主	Rabbit
特异性	Customer feedback indicates that there might be some batch-to-batch variation with respect to mouse cross-reactivity. We recommend ab124432 as an alternative antibody for mouse samples. Please contact Abcam if you have any queries.
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Pig 
免疫原	Synthetic peptide within Mouse CD31 aa 700 to the C-terminus (C terminal). The exact sequence is proprietary. Database link: Q08481
常规说明	<p>This product is FOR RESEARCH USE ONLY. For commercial use, please contact partnerships@abcam.com.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.1% Sodium azide Constituents: PBS, 1% BSA
纯度	Immunogen affinity purified

克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab28364于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P	★★★★★ (89)	<p>1/50. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.</p> <p>The ideal fixation time will depend on the size of the tissue block and the type of tissue, but fixation between 18–24h is suitable for most samples.</p> <p>Positive Control: Hu tonsil tissue</p>

靶标

功能

Induces susceptibility to atherosclerosis (By similarity). Cell adhesion molecule which is required for leukocyte transendothelial migration (TEM) under most inflammatory conditions. Tyr-690 plays a critical role in TEM and is required for efficient trafficking of PECAM1 to and from the lateral border recycling compartment (LBRC) and is also essential for the LBRC membrane to be targeted around migrating leukocytes. Prevents phagocyte ingestion of closely apposed viable cells by transmitting 'detachment' signals, and changes function on apoptosis, promoting tethering of dying cells to phagocytes (the encounter of a viable cell with a phagocyte via the homophilic interaction of PECAM1 on both cell surfaces leads to the viable cell's active repulsion from the phagocyte. During apoptosis, the inside-out signaling of PECAM1 is somehow disabled so that the apoptotic cell does not actively reject the phagocyte anymore. The lack of this repulsion signal together with the interaction of the eat-me signals and their respective receptors causes the attachment of the apoptotic cell to the phagocyte, thus triggering the process of engulfment). Isoform Delta15 is unable to protect against apoptosis. Modulates BDKRB2 activation. Regulates bradykinin- and hyperosmotic shock-induced ERK1/2 activation in human umbilical cord vein cells (HUVEC).

组织特异性

Expressed on platelets and leukocytes and is primarily concentrated at the borders between endothelial cells. Isoform Long predominates in all tissues examined. Isoform Delta12 is detected only in trachea. Isoform Delta14-15 is only detected in lung. Isoform Delta14 is detected in all tissues examined with the strongest expression in heart. Isoform Delta15 is expressed in brain, testis, ovary, cell surface of platelets, human umbilical vein endothelial cells (HUVECs), Jurkat T-cell leukemia, human erythroleukemia (HEL) and U937 histiocytic lymphoma cell lines (at protein level).

序列相似性

Contains 6 Ig-like C2-type (immunoglobulin-like) domains.

结构域

The Ig-like C2-type domains 2 and 3 contribute to formation of the complex with BDKRB2 and in regulation of its activity.

翻译后修饰


Phosphorylated on Ser and Tyr residues after cellular activation. In endothelial cells Fyn mediates mechanical-force (stretch or pull) induced tyrosine phosphorylation.

细胞定位

Membrane. Cell junction. Localizes to the lateral border recycling compartment (LBRC) and

recycles from the LBRC to the junction in resting endothelial cells and Cell junction. Localizes to the lateral border recycling compartment (LBRC) and recycles from the LBRC to the junction in resting endothelial cells.


图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD31 antibody (ab28364)

This image is courtesy of an Abreview submitted by Mr. Victor Bernal.


Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human tonsil tissue labelling CD31 with ab28364 at 1/50 dilution (30 minutes at 23°C). 10% Neutral buffered formalin used as fixative. Heat mediated antigen retrieval was performed at 95°C for 20 minutes. Blocking performed using 2.5% serum for 20 minutes at 23°C. HRP-conjugated Goat Anti-rabbit was used as the secondary antibody.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD31 antibody (ab28364)

This image is courtesy of an anonymous Abreview.

Immunohistochemistry (PFA-fixed paraffin-embedded sections) analysis of human palatine tonsil tissue labelling CD31 with ab28364 at 1/50 dilution (15 hours at 4°C). Heat mediated antigen retrieval was performed using Tris/EDTA buffer. Blocking performed with 5% normal goat serum for 1 hour at 21°C. HRP-conjugated Goat Anti-Rabbit was used as the secondary antibody.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD31 antibody (ab28364)

ab28364 at 1/50 dilution staining CD31 in human tonsil tissue section by Immunohistochemistry (Formalin/ PFA fixed paraffin embedded tissue sections).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors