

Anti-CD42d/GP5 antibody [1C2] ab95780

2 References [1 图像](#)

概述

产品名称	Anti-CD42d/GP5抗体[1C2]
描述	亚美尼亚仓鼠单克隆抗体[1C2] to CD42d/GP5
宿主	Armenian hamster
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse CD42d/GP5. (C57BL/6 x DBA/2)F1 mouse platelets
阳性对照	Mouse platelets
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: 2.68% PBS</p>
纯度	Protein G purified
克隆	单克隆
克隆编号	1C2
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab95780于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 0.125-0.25µg for 10^{5-8} cells. Staining the cell sample in a final volume of 100 µl is recommended.

靶标

功能	The GPIb-V-IX complex functions as the vWF receptor and mediates vWF-dependent platelet adhesion to blood vessels. The adhesion of platelets to injured vascular surfaces in the arterial circulation is a critical initiating event in hemostasis.
组织特异性	Platelets and megakaryocytes.
序列相似性	Contains 13 LRR (leucine-rich) repeats. Contains 1 LRRCT domain. Contains 1 LRRNT domain.
翻译后修饰	The N-terminus is blocked.
细胞定位	Membrane.

图片



Flow cytometric staining of mouse platelets with Anti-Mouse CD41 FITC and 0.125µg of Armenian Hamster IgG Isotype Control Purified (left) or 0.125µg of ab95780 (right) followed by Anti-Armenian Hamster IgG PE. Total events were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors