

Anti-CD48 antibody [MEM-102] ab9185

[2 References](#) [2 图像](#)

概述

产品名称	Anti-CD48抗体[MEM-102]
描述	小鼠单克隆抗体[MEM-102] to CD48
宿主	Mouse
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human 预测可用于: Non human primates 
免疫原	Tissue, cells or virus corresponding to Human CD48. Burkitt's lymphoma cell line Raji
阳性对照	Flow Cyt: Human peripheral blood.
常规说明	<p>This antibody is more than 95 % pure (by PAGE).</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.097% Sodium azide Constituent: PBS
纯度	Protein A purified
纯化说明	Purified from TCS. Purity >95% by SDS-PAGE.
克隆	单克隆
克隆编号	MEM-102
骨髓瘤	unknown

同种型 IgG1
轻链类型 unknown

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab9185于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 2 µg/ml.


靶标

功能 Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.


序列相似性 Contains 1 Ig-like C2-type (immunoglobulin-like) domain.
Contains 1 Ig-like V-type (immunoglobulin-like) domain.

细胞定位 Cell membrane.

图片



Flow cytometry analysis of human peripheral whole blood labeling CD48 with ab9185 at 3 µg/mL (GAM APC).



Flow cytometry analysis of human peripheral whole blood labeling CD48 positive lymphocytes with ab9185 at 3 µg/mL (GAM APC) from neutrophil granulocytes (black-dashed).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors