

Anti-CD5 antibody [CC17] ab34070

[1 References](#) [1 图像](#)

概述

产品名称	Anti-CD5抗体[CC17]
描述	小鼠单克隆抗体[CC17] to CD5
宿主	Mouse
特异性	ab34070 is specific for bovine CD5. This clone reacts with BoCD5.1, which is a polymorphic antigen expressed on cells of all Bos Taurus animals and a small proportion of Bos Indicus animals.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Cow
免疫原	Tissue, cells or virus corresponding to Cow CD5. Bovine thymocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide
纯度	Protein G purified
克隆	单克隆
克隆编号	CC17
骨髓瘤	NS1
同种型	IgG1

The Abpromise guarantee **Abpromise™**承诺保证使用ab34070于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

应用说明

FACS: 1/25 - 1/100. Use 10µl of the working dilution for 10⁶ cells in 100ul.
 IHC-Fr: Use at an assay dependent dilution.
 IP: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标**功能**


May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.

序列相似性

Contains 3 SRCR domains.

细胞定位

Membrane.

图片

Staining of bovine peripheral blood lymphocytes with Mouse anti Bovine CD5.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors