

Anti-CD57 antibody [TB01] ab33273

1 References

概述

产品名称	Anti-CD57抗体[TB01]
描述	小鼠单克隆抗体[TB01] to CD57
宿主	Mouse
经测试应用	适用于: IHC-Fr, IHC-P, Flow Cyt 不适用于: IP or WB
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD57. Human neuroblastoma cells
阳性对照	Human tonsil tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Ion Exchange Chromatography
克隆	单克隆
克隆编号	TB01
骨髓瘤	P3x63-Ag8.653
同种型	IgM

应用

The Abpromise guarantee Abpromise™承诺保证使用ab33273于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		1/10 - 1/100.
IHC-P		1/10 - 1/50. This product does not require protein digestion pre-treatment of paraffin sections. This product does not require antigen retrieval using heat treatment prior to staining of paraffin sections but sodium citrate buffer may enhance staining.
Flow Cyt		1/50 - 1/100. ab91545 - Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody. Use 10µl of the suggested working dilution to label 10 ⁶ cells in 100µl.

应用说明 Is unsuitable for IP or WB.

靶标

相关性 CD57 is a carbohydrate epitope that contains a sulfoglucuronyl residue in several adhesion molecules of the nervous system. Is created by Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 (B3GAT1) The CD57 antigen is present on subsets of T cells and natural killer cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors