

Anti-CD6/T12 antibody [MEM-98] ab9181

2 References

概述

产品名称	Anti-CD6/T12抗体[MEM-98]
描述	小鼠单克隆抗体[MEM-98] to CD6/T12
宿主	Mouse
特异性	CD6 / T12 antigen (human)
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Full length native protein (purified) corresponding to Human CD6/T12.
阳性对照	Spleen.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Protein A purified
纯化说明	Purified from TCS. Purity >95% by SDS-PAGE.
克隆	单克隆
克隆编号	MEM-98
骨髓瘤	unknown
同种型	IgG1

轻链类型

unknown

应用

The Abpromise guarantee Abpromise™承诺保证使用ab9181于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 2 µg/ml.

靶标

功能	Involved in cell adhesion. Binds to CD166.
组织特异性	Expressed by thymocytes, mature T-cells, a subset of B-cells known as B-1 cells, and by some cells in the brain.
序列相似性	Contains 3 SRCR domains.
翻译后修饰	After T-cell activation, becomes hyperphosphorylated on Ser and Thr residues and phosphorylated on Tyr residues. Contains intrachain disulfide bond(s).
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors