

Anti-CD73 antibody [7G2] ab54217

★★★★★ [3 Abreviews](#) [29 References](#) [2 图像](#)


概述

| | |
|-------|---|
| 产品名称 | Anti-CD73抗体[7G2] |
| 描述 | 小鼠单克隆抗体[7G2] to CD73 |
| 宿主 | Mouse |
| 经测试应用 | 适用于: Flow Cyt, ICC 不适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Full length purified Human placental CD73. |
| 阳性对照 | Flow Cyt: Human peripheral blood lymphocytes. ICC: U-87 MG cells. |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| 存储溶液 | pH: 7.40 Preservative: 0.1% Sodium azide Constituent: PBS |
| 纯度 | Protein A purified |
| 克隆 | 单克隆 |
| 克隆编号 | 7G2 |
| 同种型 | IgG2a |
| 轻链类型 | kappa |

| 应用 | Ab评论 | 说明 |
|----------|-----------|--|
| Flow Cyt | ★☆☆☆☆ (1) | Use a concentration of 1 - 2 µg/ml. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody. |
| ICC | ★★★★★ (1) | 1/100. |



Flow Cytometry - Anti-CD73 antibody [7G2]
(ab54217)

ab54217 staining human peripheral blood lymphocytes by FACS analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors